




GYRO ST-ER : GYRO ER collet center alignment device with a cylindrical shank.



Designation	SS	CSI	d min	d max	D3	D1	D2	L1	L	L2	Kg
GYRO ST20 ER20	20	ER20	1.0	13.0	34.00	57.0	63.00	28.5	58.80	80.00	0.84
GYRO ST20 ER25	20	ER25	1.0	16.0	34.00	57.0	63.00	28.5	58.80	80.00	1.40
GYRO ST25 ER25	25	ER25	1.0	16.0	42.00	74.0	79.00	35.5	65.65	80.00	1.46
GYRO ST25 ER32	25	ER32	2.0	20.0	50.00	74.0	79.00	36.5	66.65	80.00	1.51
GYRO ST32 ER32	32	ER32	2.0	20.0	50.00	74.0	79.00	36.5	66.65	80.00	1.70
GYRO ST40 ER32	40	ER32	2.0	26.0	50.00	74.0	79.00	36.5	66.65	80.00	2.02

For adaptation Click on "More Info".

Radial adjustment 2.0 mm. Angular adjustment 1°.

First-time users should buy a GYRO kit which includes a test bar and a bushing for performing the alignment procedure.